## Evaluation and recommendation of sample preservation techniques for 16S rRNA sequencing in oilfield systems

Jordan J. Schmidt<sup>1</sup>, Brad McIlwain<sup>1</sup>, Wei Huang<sup>1</sup>, Gilberto Guerrero<sup>2</sup>, Prasad Dhulipala<sup>2</sup>, Mary Eid<sup>2</sup>, Neil Sharma<sup>1</sup> <sup>1</sup>LuminUltra Technologies Ltd, 520 King Street, Fredericton, NB, Canada, E3B 6G3 <sup>2</sup>Baker Hughes, Houston, TX

# LUMINUTRA microbial monitoring



#### BACKGROUND

The implementation of metagenomics in oil and gas industry has drastically increased in the last several years. It is being used for a host of applications, including: detection of the origin of microbially influenced corrosion (MIC), biocide efficacy studies, and microbial content evaluations of injection water and water used for hydraulic fracturing. One major deficiency with the application of metagenomics, is the lack of standardized protocols. Due to the remote nature of oilfields, sample preservation and maintaining integrity of microbial population during transport is of particular concern. If not preserved properly, the microbial community composition can change between sample collection and DNA extraction.



#### OBJECTIVE

To evaluate several common preservation techniques that are suitable for remote oilfield locations and contribute towards the development of standard protocols for metagenomics analysis.

#### METHODOLOGY

Two samples were collected from an oilfield and shipped to the lab, the DNA was extracted from the control sample. Samples were then preserved using various methods and subjected to different scenarios including a simulated shipping delay and temperature cycling. Preserved samples were held for 7 and 14 days prior to DNA extraction and analysis. All samples were analyzed for total prokaryotes using qPCR and sequenced using Illumina 16s rRNA sequencing (515/806 primer).

	<section-header><section-header></section-header></section-header>	<ol> <li>Control</li> <li>LuminUltra GeneCount Chemical Preservation (7, 14 days)</li> <li>Filtered and dried (7, 14 days)</li> </ol>	qPCR – Total 16S rRNA Prokaryotes Sequencing
		<ul> <li>4. Unpreserved (7, 14 days) Temperature cycled (12 h at 20°C, 12h at 40°C)</li> <li>5. Refrigerated (7, 14 days)</li> <li>6. Frozen (14 days)</li> <li>7. Shipping delay (Refrigerated for 24 hours, room temperature for 3 days)</li> </ul>	
	Sample Collection	Sample Preservation	Sample Analysis

Two samples were also preserved in the field using LuminUltra's GeneCount chemical preserved samples were shipped to the lab where they were further temperature cycled and analyzed after 7 and 14 days.

#### RESULTS – IN LAB

Well Water

#### RESULTS-IN FIELD

Well Water



#### OUTCOMES – IN LAB

- There were three preservation techniques that gave similar results to the control in both samples: freezing, refrigeration and GeneCount chemical preservation.
- Filtering and drying had similar results to the control in well water (lower alpha diversity) but not in frac pit water (higher alpha diversity), which shows that it does not universally preserve all samples. There was a considerable increase in Bacilli-associated genera in the filtered frac pit water.
- The GeneCount chemical preservation kit did result in some additional low abundance OTUs, 3. especially in the highly turbid well water. It is hypothesized that this is due to the extended extraction period offered by chemical preservation which resulted in the lysing of additional and more robust cells.

Low Sim	ilarity	High Similarity	
Filtered-14d	Filtered-7d	GeneCount-14d	GeneCount-7d

- Well Water Similar to the laboratory results, filtration and GeneCount chemical preservation gave similar results for both hold times. Compared to each other the preservation methods
- Frac Pit Water The GeneCount chemical preservation kit gave similar results after 7 and 14 days of temperature cycling. Like the laboratory results, the filtered samples had low similarity when compared to each other and compared to the chemically preserved samples. The filtered samples were characterized by relative increases in Pseudomonas (day 7), Bacillus (day 7 and 14) and Streptomyces (day 7 and 14) compared to the chemically preserved samples.

### CONCLUSIONS

Various preservation methods were tested on two oilfield sample matrices. Three preservation methods performed well for both samples: freezing, refrigeration and GeneCount chemical preservation. Filtration and air drying worked well for one sample, but not the other. This shows that testing multiple sample types is important when evaluating preservation methods.

The results show that chemical sample preservation is recommended, particularly in cases where cold-chain transport is not possible.