

PRODUCT FACT SHEET

DIS—LBUN and DIS—SBUN

Disposable Equipment Bundles



Product Overview

As with many common laboratory measurements, LuminUltra's test kits require the use of reliable liquid transfer equipment and mixing vessels. For flexibility and convenience, LuminUltra recommends the use of disposable test tubes and micropipette tips for analyzing samples with its test kits to minimize cost, maintenance, and the potential of analysis contamination. Both **DIS-LBUN** and **DIS-SBUN** meet this requirement and are applicable to all test kit models, thus requiring small investment with great flexibility!



Product Compatibility

For your convenience, this product is compatible with all 5 LuminUltra test kits.

Catalogue Number	Designed for...
DIS-SBUN	All 12 or 25 size test kits
DIS-LBUN	All 50/100 or 125/250 size test kits

Product Features

This bundle of equipment has the following features:

Included in DIS-SBUN/LBUN: (Individual Cat. #)	Quantity		Compatible with: (Individual Cat. #)
	SBUN	LBUN	
17x100mm Polypropylene Culture Tubes (Cat. #: DIS-CT17)	50	500	17mm Culture Tube Rack (Cat #: EQP-TR17)
17mm Polypropylene Culture Tube Caps (Cat.# DIS-C17)	50	500	17x100mm Polypropylene Culture Tubes (Cat #: DIS-CT17)
12x55mm Polypropylene Culture Tubes (Cat. #: DIS-CT12)	50	500	12mm Culture Tube Rack (Cat #: EQP-TR12) Kikkoman C-100 LumiTester (Cat #: EQP-LUK)
100-1000µL (Blue) Polypropylene Pipette Tips (Cat. #: DIS-PT1)	100	1000	100-1000µL Adjustable Micropipettor (Cat #: EQP-MP1)
1-5mL (Green) Polypropylene Pipette Tips (Cat. #: DIS-PT5)	40	400	1-5mL Adjustable Micropipettor (Cat #: EQP-MP5)

*The disposable equipment components in **DIS-SBUN** and **DIS-LBUN** are designed to be specifically compatible with the Kikkoman C-100 LumiTester (**EQP-LUK**) and our Laboratory Equipment Bundle (**EQP-BUN**). If you require a different configuration, please consult LuminUltra.



ISO 9001:2000

LuminUltra™ Technologies Ltd.
Innovative Biological Monitoring Technologies

440 King Street, King Tower, Suite 630, Fredericton, New Brunswick, Canada E3B 5H8

Tel: (506) 459-8777 Fax: (506) 453-9860 www.luminultra.com